

**Mianserin-D3 solution, 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-M-901-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|--------------|---------------------|-------|--------|--------|
|                  | Mianserin-D3      | [81957-76-8] | 100                 | µg/mL |        |        |