

Dehydroepiandrosterone-D5 (DHEA-D5) (2,2,3,4,4-D5) solution, 100 µg/mL in methanol, ampule of 1 mL, certified reference material

Art. ID SAF-D-064-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium Prasterone Sulfa te-d4	[97453-25-3]	100	µg/mL		