

**Dehydroepiandrosterone-D5-3-sulfate (DHEAS-D5) (2,2,3,4,4,-D5) sodium salt solution, 100 µg/mL in methanol (as free sulfate), ampule of 1 mL, certified reference material**

Art. ID                      SAF-D-066-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dehydroepiandrosterone- D5-3-sulfate (DHEAS-D5) (2,2,3,4,4,-D5) sodium salt		100	µg/mL		