

**N-Desmethyldoxepin-D3 hydrochloride solution, cis/trans, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material**

Art. ID                      SAF-D-075-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Desmethyldoxepin-D3 H Cl(cis/trans)		100	µg/mL		