

**(±)-4-Hydroxy-3-methoxymandelic Acid-D3 (ring-D3) solution, 100 µg/mL in methanol, ampule of 1 mL, certified reference material**

Art. ID                      SAF-H-091-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(±)-4-Hydroxy-3-methoxy mandelic acid-D3 (ring- D3)		100	µg/mL		