

(+/-)-Norpseudoephedrine-D3 hydrochloride solution, 100 µg/mL in methanol (as free base), ampule of 1 mL, certified reference material

Art. ID SAF-N-087-1ML
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(±)-Norpseudoephedrine-D3 hydrochloride		100	µg/mL		
	(±)-Norpseudoephedrine-D3 hydrochloride		100	µg/mL		