

**Anti-ARHGEF1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA012924-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-115 kDa guanine nucleotide exchange factor|Anti-Rho guanine nucleotide exchange factor  
1|Anti-Sub15|Anti-p115-RhoGEF|Anti-p115RhoGEF|