

**Anti-PCOLCE2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA013203-25UL  
Unit                         1 x 25 µL  
Deliverydetails           /not restricted

Description

Synonyms: |Anti-PCPE-2|Anti-Procollagen C-endopeptidase enhancer 2 precursor|Anti-Procollagen C-proteinase enhancer 2|Anti-Procollagen COOH- terminal proteinase enhancer 2|