

Anti-GGT7 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA013204-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-GGT 4|Anti-Gamma-glutamyltransferase 4 precursor|Anti-Gamma-glutamyltransferase-like 3|Anti-Gamma-glutamyltransferase-like 5|Anti-Gamma-glutamyltranspeptidase 4|