

Anti-FCRL1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA013323-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Fc receptor homolog 1|Anti-Fc receptor-like protein 1 precursor|Anti-FcR-like protein 1|Anti-FcRH1|Anti-FcRL1|Anti-IFGP family protein 1|Anti-Immune receptor translocation-associated protein 5|Anti-hIFGP1|