

Anti-FKBP14 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA013329-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-22 kDa FK506-binding protein|Anti-FK506-binding protein 14 precursor|Anti-FKBP-22|Anti-PPIase|Anti-Peptidyl-prolyl cis-trans isomerase|Anti-Rotamase|