

Anti-MYL7 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA013331-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-MLC-2a|Anti-MLC2a|Anti-Myosin light chain 2a|Anti-Myosin regulatory light chain 2, atrial isoform|Anti-Myosin regulatory light chain 7|