

Anti-SLC13A3 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA014736-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Na(+)/dicarboxylate cotransporter 3|Anti-NaDC-3|Anti-Sodium-dependent high-affinity dicarboxylate transporter 2|Anti-Solute carrier family 13 member 3|Anti-hNaDC3|