

**Anti-TIMD4 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA015625-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-T-cell immunoglobulin and mucin domain-containing protein 4 precursor|Anti-T-cell membrane protein 4|Anti-TIM-4|Anti-TIMD-4|