

Anti-PSMB6 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023312-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Macropain delta chain|Anti-Multicatalytic endopeptidase complex delta chain|Anti-Proteasome delta chain|Anti-Proteasome subunit Y|Anti-Proteasome subunit beta type-6|