

Anti-BRF2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023378-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-B-related factor 2|Anti-BRF-2|Anti-TFIIIB50|Anti-Transcription factor IIIB 50 kDa subunit|Anti-hBRFU|Anti-hTFIIIB50|