

Anti-USP43 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023389-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Deubiquitinating enzyme 43|Anti-Ubiquitin carboxyl-terminal hydrolase 43|Anti-Ubiquitin thioesterase 43|Anti-Ubiquitin-specific-processing protease 43|