

Anti-AGAP2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023474-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-AGAP-2|Anti-Arf-GAP, GTPase, ANK repeat and PH domain-containing protein
2|Anti-Centaurin-gamma-1|Anti-GGAP2|Anti-GTP-binding and GTPase-activating protein
2|Anti-PIKE|Anti-Phosphatidylinositol-3-kinase enhancer|