

**Anti-ADPRM antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution, Ab2**

Art. ID SAF-HPA023820-25UL  
Unit 1 x 25 µL  
Deliverydetails /not restricted

Description

Synonyms: |Anti-ADPRibase-Mn|Anti-C17orf48|Anti-CDP-choline phosphohydrolase|Anti-Magnesium-dependent ADP-ribose/CDP-alcohol diphosphatase|