

Anti-DMD antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023885-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms:

|Anti-BMD|Anti-DXS142|Anti-DXS164|Anti-DXS206|Anti-DXS230|Anti-DXS239|Anti-DXS268|Anti-DXS269|Anti-DXS
270|Anti-DXS272|Anti-MRX85|