

Anti-ARHGAP27 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA023919-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-ARHGAP27|Anti-CIN85-associated multi-domain-containing Rho GTPase-activating protein 1|Anti-Rho GTPase-activating protein 27|Anti-Rho-type GTPase-activating protein 27|Anti-SH3D20|