

Anti-ATP6V0D1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA016938-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-32 kDa accessory protein|Anti-P39|Anti-V-ATPase 40 kDa accessory protein|Anti-V-ATPase AC39 subunit|Anti-V-ATPase subunit d 1|Anti-Vacuolar proton pump subunit d 1|