

Anti-EMCN antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID	SAF-HPA005928-25UL
Unit	1 x 25 µL
Deliverydetails	/not restricted

Description

Synonyms: |Anti-Endomucin precursor antibody produced in rabbit|Anti-Endomucin-2 antibody produced in rabbit|Anti-Gastric cancer antigen Ga34 antibody produced in rabbit|Anti-Mucin-14 antibody produced in rabbit|