

Anti-WTAP antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010550-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Female-lethal(2)D homolog|Anti-Pre-mRNA-splicing regulator WTAP|Anti-WT1-associated protein|Anti-Wilms tumor 1-associating protein|Anti-hFL(2)D|