

Anti-NAXD antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010551-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-CARKD|Anti-FLJ10769|Anti-LP3298|Anti-RP11-90L1.3 antibody produced in rabbit|Anti-YM014_HUMAN Isoform 2 of Q8IW45 - Homo sapiens (human)|