

Anti-IL6ST antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010558-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-CD130 antigen antibody produced in rabbit|Anti-CDw130 antibody produced in rabbit|Anti-IL-6R-beta antibody produced in rabbit|Anti-Interleukin-6 receptor subunit beta precursor antibody produced in rabbit|Anti-Interleukin-6 signal transducer anti|