

Anti-NT5C3A antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010630-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms:

|Anti-NT5C3|Anti-P5'N-1|Anti-PN-I|Anti-POMP|Anti-PSN1|Anti-UMPH|Anti-UMPH1|Anti-cN-III|Anti-hUMP1|Anti-p36|