

**Anti-PIANP antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA010631-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-C12orf53|Anti-DKFZp547D2210|Anti-H3/j|Anti-H3FJ|Anti-Uncharacterized protein C12orf53  
Precursor|