

**Anti-PTGER3 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA010689-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-PGE receptor, EP3 subtype|Anti-PGE2-R|Anti-Prostaglandin E2 receptor EP3 subtype|Anti-Prostanoid EP3 receptor|