

**Anti-ALK antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA010694-25UL  
Unit                         1 x 25 µL  
Deliverydetails           /not restricted

Description

Synonyms: |Anti-ALK tyrosine kinase receptor precursor antibody produced in rabbit|Anti-Anaplastic lymphoma kinase antibody produced in rabbit|Anti-CD246 antigen antibody produced in rabbit|