

Anti-PRRG2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010702-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-Proline-rich gamma-carboxyglutamic acid protein 2 antibody produced in rabbit|Anti-Proline-rich Gla protein 2 antibody produced in rabbit|Anti-Transmembrane gamma-carboxyglutamic acid protein 2 precursor antibody produced in rabbit|