

Anti-FCGR2A antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010718-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-CD32 antigen antibody produced in rabbit|Anti-CDw32 antibody produced in rabbit|Anti-Fc-gamma RII-a antibody produced in rabbit|Anti-Fc-gamma-RIIa antibody produced in rabbit|Anti-FcRII-a antibody produced in rabbit|Anti-IgG Fc receptor II-a ant|