

Anti-C1orf43 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010725-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-HCV NS5A-transactivated protein 4 antibody produced in rabbit|Anti-Hepatitis C virus NS5A-transactivated protein 4 antibody produced in rabbit|Anti-Protein NICE-3 antibody produced in rabbit|Anti-S863-3 antibody produced in rabbit|Anti-Uncharacte|