

Anti-DCT antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010743-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-DCT|Anti-DT|Anti-L- dopachrome Delta-isomerase|Anti-L-dopachrome tautomerase precursor|Anti-TRP-2|Anti-TRP2|Anti-Tyrosinase-related protein 2|