

**Anti-B4GALT3 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID SAF-HPA010793-25UL  
Unit 1 x 25 µL  
Deliverydetails /not restricted

Description

Synonyms: |Anti-beta-1,4-GalTase 3 antibody produced in rabbit|Anti-beta-1,4-Galactosyltransferase 3 antibody produced in rabbit|Anti-beta4Gal-T3 antibody produced in rabbit|Anti-UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 3 antibody produced in rabbit|An|