

Anti-B4GALT1 antibody produced in rabbit, Ab1, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010806-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-beta-1,4-GalTase 1 antibody produced in rabbit|Anti-beta-1,4-Galactosyltransferase 1 antibody produced in rabbit|Anti-beta4Gal-T1 antibody produced in rabbit|Anti-UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 1 antibody produced in rabbit|An|