

Anti-DNAJB11 antibody produced in rabbit, Ab1, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA010814-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-DnaJ homolog subfamily B member 11 precursor antibody produced in rabbit|Anti-ER-associated Hsp40 co-chaperone antibody produced in rabbit|Anti-ER-associated dnaJ protein 3 antibody produced in rabbit|Anti-ErJ3 antibody produced in rabbit|Anti-PW|