

**Copper base, Type: Bronze - Aluminum, RM, Mushroom 60 mm x 5 mm**

Art. ID CTIF-2151-R  
Unit each  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 9,43                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 84,75               | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 4,48                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,73                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,56                | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | <0,005              | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,015               | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | <0,01               | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | <0,01               | %    |        |        |