

**Anti-SLC11A1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA029590-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-Lsh|Anti-Nramp|Anti-Nramp1|Anti-Solute carrier family 11 (proton-coupled divalent metal ion transporters), member 1|