

**Anti-ACAP2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA034807-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-ArfGAP with coiled-coil, ankyrin repeat and PH domains  
2|Anti-CENTB2|Anti-CNT-B2|Anti-KIAA0041|