

**Anti-PDCD4 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA027214-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-Neoplastic transformation inhibitor protein|Anti-Nuclear antigen H731-like|Anti-Programmed cell death protein 4|Anti-Protein 197/15a|