

**Anti-MAGED2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution, ab1**

Art. ID                      SAF-HPA031572-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms:

|Anti-11B6|Anti-BCG1|Anti-HCA10|Anti-JCL-1|Anti-MAGE-D2|Anti-MAGED|Anti-MGC8386|Anti-melanoma antigen family D, 2|