

Anti-STRADA antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA031637-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-LYK5|Anti-MGC14832|Anti-NY-BR-96|Anti-ORB3|Anti-STE20-related kinase adaptor alpha|Anti-STRAD|Anti-Stlk|Anti-XTP4|