

Anti-ZFP64 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA035112-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-DJ548G19.1|Anti-DJ831D17.1|Anti-FLJ10734|Anti-FLJ12628|Anti-ZNF338|Anti-Zinc finger protein 64 homolog (mouse)|