

**Anti-ATP2C1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID SAF-HPA035116-25UL  
Unit 1 x 25 µL  
Deliverydetails /not restricted

Description

Synonyms: |Atp2C1 Antibody|Atp2C1 Antibody - Anti-ATP2C1 antibody produced in rabbit|Anti-ATP2C1A|Anti-ATPase, Ca<sup>++</sup> transporting, type 2C, member 1|Anti-BCPM|Anti-KIAA1347|Anti-PMR1|Anti-SPCA1|