

**Copper base, Type: Bronze - Aluminum, RM, Mushroom 60 mm x 5 mm**

Art. ID CTIF-3610-Q  
Unit each  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 7,1                 | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,09                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 82,32               | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 3,98                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,045               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 5,4                 | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,23                | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,065               | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,25                | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 0,51                | %    |        |        |