

Copper base, Type: Bronze - Aluminum, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-4149-G
Unit each
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 9,84 | % | | |
| | Copper (Cu) | [7440-50-8] | 84,95 | % | | |
| | Iron (Fe) | [7439-89-6] | 2 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,21 | % | | |
| | Nickel (Ni) | [7440-02-0] | 1,96 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,15 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,18 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,34 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,37 | % | | |