

Copper base, Type: Beryllium-Cobalt, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-4584
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,033	%		
	Beryllium (Be)	[7440-41-7]	2,53	%		
	Cobalt (Co)	[7440-48-4]	0,04	%		
	Copper (Cu)	[7440-50-8]	97,05	%		
	Iron (Fe)	[7439-89-6]	0,12	%		
	Nickel (Ni)	[7440-02-0]	0,015	%		
	Silicon (Si)	[7440-21-3]	0,166	%		
	Tin (Sn)	[7440-31-5]	0,022	%		
	Zinc (Zn)	[7440-66-6]	0,022	%		