

**Copper base, Type: Beryllium-Cobalt, RM, Mushroom 60 mm x 5 mm**

Art. ID                      CTIF-4873  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,094	%		
	Beryllium (Be)	[7440-41-7]	0,17	%		
	Cobalt (Co)	[7440-48-4]	0,98	%		
	Chromium (Cr)	[7440-47-3]	0,105	%		
	Copper (Cu)	[7440-50-8]	98,4	%		
	Iron (Fe)	[7439-89-6]	0,078	%		
	Nickel (Ni)	[7440-02-0]	0,049	%		
	Silicon (Si)	[7440-21-3]	0,088	%		
	Tin (Sn)	[7440-31-5]	~0,007	%		
	Zinc (Zn)	[7440-66-6]	~0,003	%		