

Copper base, Type: Bronze, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-B10
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,205	%		
	Arsenic (As)	[7440-38-2]	0,0075	%		
	Copper (Cu)	[7440-50-8]	83,65	%		
	Iron (Fe)	[7439-89-6]	0,165	%		
	Nickel (Ni)	[7440-02-0]	1,01	%		
	Phosphorus (P)	[7723-14-0]	0,014	%		
	Lead (Pb)	[7439-92-1]	4,07	%		
	Sulfur (S)	[7704-34-9]	0,05	%		
	Antimony (Sb)	[7440-36-0]	1,14	%		
	Tin (Sn)	[7440-31-5]	6,95	%		
	Zinc (Zn)	[7440-66-6]	2,75	%		